

MODEL/PART NUMBER	PT-2	PT-10	PT-30	PT-60	NEW! PTR 7.5	NEW! PTR 75
Transmitter Type	4-wire					
PendoTECH Single-Use Pressure Sensor	Yes					
Number of Sensor Channels	1 pressure channel					
Installation & Enclosure Dimension	DIN Rail Mount; Request drawing					
Power Input	Nominal 24VDC; 9-36V allowed					
Pressure Input	0-2psi (0.138bar)	10psi (.69bar)	30psi (2.07bar)	60psi (4.14bar)	-11.5 - 7.5psi	-11.5 - 75psi
Reusable Sensor Cable Installed	Yes				Cable in package; Doesn't come pre-installed	
Relays	-				(1) NO/NC	
Outputs/Communication	(1) Active 4-20mA output				(1) Active 4-20mA output & Modbus-TCP	
Approvals	CE Mark EN61000-6-2:2015; EN61000-6-4:2011; EN61010-1:2010; EN13849-1:2015; REACH and RoHS Compliant; UL and cUL Recognized				CE & UKCA Mark EN IEC 61326-1:2021; EN 61010-1:2010/A1:2019; EN ISO 13849-1:2023 RoHS and REACH compliant; UL and cUL Recognized	