

MODEL/PART NUMBER	PTR 7.5	PTR 75
Transmitter Type	4-wire	
Pendotech Single-Use Pressure Sensor	Yes	
Number of Sensor Channels	1 pressure channel	
Installation & Enclosure Dimension	DIN Rail Mount; Request drawing	
Power Input	Nominal 24VDC; 9-36V allowed	
Pressure Input	-11.5 - 7.5psi	-11.5 - 75psi
Reusable Sensor Cable Installed	Cable in package; Doesn't come pre-installed	
Relays	(1) NO/NC	
Outputs/Communication	(1) Active 4-20mA output & Modbus-TCP	
Approvals	CE & UKCA Mark EN IEC 61326-1:2021; EN 61010-1:2010/A1:2019; EN ISO 13849-1:2023 RoHS and REACH compliant; UL and cUL Recognized	